



Seat Colors  
Black Vinyl



Body Colors  
Red  
Blue



**HD Version  
450 lb.  
weight capacity!**

golden  
**COMPASS™**  
**HD**  
GP620

The heavy-duty version of the original Compass offers users a robust 450 lb. weight capacity with all of the same great features including the patented springless articulating chassis and true center-wheel drive design for unmatched maneuverability in tight spaces. The Compass HD comes in a variety of seat sizes up to 24" wide.



## Golden Compass™ HD Specifications

<b>Model</b>	<b>GP620</b>
Model Name	Compass HD
Medicare Code	K0825
Weight Capacity	450 lbs
Drive Wheels	Center
Maximum Speed	4 mph
Operating Range <sup>1</sup>	12 Miles <sup>3</sup>
Ground Clearance	3.5"
Turning Radius	20"
Type Batteries	2-NF22
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	75 Amp
<b>Dimensions:</b>	
Length	40"
Width	25"
Height (ground to top of back of seat)	38.25" - 42.75"
Height (ground to top of back of headrest)	47.25" - 53.75"
Ground to Top of Seat	21.25" - 24.25"
Number of Seat Height Positions	5
Size of Increments	0.5"
<b>Weight of Unit: (assembled)</b>	226 lbs
Base	101 lbs
Seat w/Arms	51 lbs
Batteries	37 lbs ea.
<b>Tires:</b>	Flat Free
Center	10"
Casters	6"
<b>Standard Seat:</b>	Captain's Seat
Back Height - no headrest	19"
Back Height with headrest	26.5" - 30"
Width x Depth	20" x 20"/18"
Color	Black Vinyl
<b>Optional Seats:</b>	Captain's Seat
Back Height - no headrest	19"
Back Height with headrest	26.5" - 30"
Width x Depth	18", 20", 22" or 24" x 18" or 20"
Color	Black Vinyl
<b>Footrest:</b>	Yes
Height adjustable	3" - 7"
Angle adjustable	0 - 35 degrees
<b>Warranty:</b>	
Frame	5 Years
Drive Train	13 Months
Electronics <sup>2</sup>	13 Months
Batteries	12 Months

<sup>1</sup> Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain & battery type.

<sup>2</sup> Electronics warranty excludes batteries.

<sup>3</sup> Battery range at maximum weight capacity and will vary due to rider weight, drive surface, terrain & battery type.



## Power Chair Accessories and Options

- Captain's seat
- Swing-away legrests
- Elevating legrests
- Stump support
- Footplate extender
- Cane/crutch holder kit
- Walker holder kit
- Oxygen tank holder kit
- Pack N Go

## 2 exciting colors



Red



Blue