

M I D S I Z E S C O O T E R



Seat Color
Black Vinyl

Body Colors
Red
Blue

GOLDEN COMPANION

3 wheel (Midsize)
GC240C

The Companion scooter provides great style with high performance and user-friendly disassembly. New features include an adjustable LED headlight; new ergonomic control panel design with backlit battery gauge display; and wraparound Delta tiller. The Companion offers great value in a full-time portable scooter.



Companion 3 wheel Specifications

Model	GC240C
Medicare Code	K0807
Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	5.3 mph
Operating Range ¹	19.9 miles
Ground Clearance Under Deck	3.75"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	37"
Type Batteries	2-U1
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	90 Amp Dynamic R Series
Dimensions:	
Length	42.25"
Width	23.5"
Height (ground to top of back of seat)	37" - 39"
Ground to Top of Seat	23" - 25"
Ground to Top of Deck	6.25"
Deck to Top of Seat	17.5" - 19.5"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	30"
Weight of Unit: (assembled)	163 lbs
Front Half	45 lbs
Rear Half	48 lbs
Seat w/Arms	26 lbs
Batteries	21 lbs ea
Tires:	Solid, Low Profile
Front	9"
Rear	10.75"
Rear Anti-Tippers	3"
Standard Seat:	Stadium Style
Slide Seat	Standard
Back Height - no headrest	16"
Width x Depth	18" x 16"
Color	Black Vinyl
Optional Seats:	Stadium Style
Back Height - no headrest	16"
Width x Depth	20" x 18"
Color	Black Vinyl
Seat Options:	
Pneumatic Rise Seat	Yes
Power Elevating Seat	Yes
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Batteries	1 Year



Angle adjustable LED headlight creates a wide path of bright light



EZ Reach tiller adjustment



Patented infinite adjustable tiller




Quiet & maintenance free motor

Standard Features on Every Model:

- Under mat troubleshooting guide
- Delta tiller with auto lock mechanism
- Easy reach tiller adjustment
- Adjustable sliding seat
- Adjustable LED headlight
- U-1 batteries standard



Disassembles easily!



Blue Red

¹. Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

². Electronics warranty excludes batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

