

FULL SIZE SCOOTER



Seat Color
Black Vinyl

Body Colors
Red
Blue

GOLDEN COMPANION

4 wheel
GC440C

The Companion four-wheel scooter offers the most leg room and foot room in its class with a modern, rugged design. Features include easy disassembly; an adjustable LED headlight; ergonomic control panel design with backlit battery gauge display; and wraparound Delta tiller. The Companion four-wheel is rugged, dependable and an outstanding value for a full-time, four-wheel scooter.



GC440C

Companion 4 wheel Specifications

Model	GC440C
Medicare Code	K0807
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	5.3 mph
Operating Range ¹	18.9 miles (U1) 23.4 miles (40AH)
Ground Clearance Under Deck	4.5"
Ground Clearance to Transaxle/Motor	2.5"
Turning Radius	56.5"
Type Batteries	2-U1, (2-40AH optional)
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	90 Amp Dynamic R Series
Dimensions:	
Length	47.5"
Width	23.5"
Height (ground to top of back of seat)	37" - 39"
Ground to Top of Seat	24" - 26"
Ground to Top of Deck	6.25"
Deck to Top of Seat	18" - 20"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	35.5"
Weight of Unit: (assembled)	
	194 lbs (U1) 210 lbs (40AH)
Front Half	74 lbs
Rear Half	48 lbs
Seat w/Arms	26 lbs
Batteries	21 lbs ea. (U1), 29 lbs ea. (40AH)
Tires: Solid/Low Profile	
Front	10.75"
Rear	10.75"
Rear Anti-Tippers	3"
Standard Seat ³: Stadium Style	
Slide Seat	Standard
Back Height - no headrest	16"
Width x Depth	18" X 16"
Color	Black Vinyl
Optional Seats: Stadium Style	
Back Height - no headrest	16"
Width x Depth	20" x 18"
Color	Black Vinyl
Optional Seats 2: Captain's Seat	
Back Height - no headrest	19"
Back Height with headrest	29"
Width x Depth	18" x 18"
Color	Black Vinyl
Seat Options:	
Pneumatic Rise Seat	Yes
Power Elevating Seat	Yes
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Batteries	1 Year



Angle adjustable LED headlight creates a wide path of bright light



EZ Reach tiller adjustment




Patented infinite adjustable tiller



Quiet & maintenance free motor

Standard Features on Every Model:

- Under mat troubleshooting guide
- Delta tiller with auto lock mechanism
- Easy reach tiller adjustment
- Adjustable sliding seat
- Adjustable LED headlight
- U-1 batteries standard



Blue Red



Disassembles easily!

¹ Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.
² Electronics warranty excludes batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.