



TRIAXE Sport



- **Rust Free Aluminum Construction**
- **Durable – 350 Weight Capacity**
 - **Ultra-Lightweight (44lbs)**
 - **Fast - Speeds Up to 12mph**
- **Up to 35 Miles on a Single Charge**
- **Folds Flat Effortlessly for Storage**
 - **Deluxe Basket Included**
- **6 Exciting Colors to Choose From**

Model	Triaxe Sport T3045
Overall Length	34"
Overall Width	20 1/2 "
Tires	8" Solid Front Tire 9.5 Solid Rear Tires
Maximum Speed	Up to 12mph
Range	Up to 35 miles
Turning Radius	33"
Weight	44lbs (without battery)
Color	Red, Light Blue, Purple, Orange, Black, Silver (limited edition colors)
Brake	Hand brake controlling rear disc brakes
Drive System	Front Wheel Drive, Brushless Motor
Handlebar	Aluminum with Rubber Hand Grips
Maximum Weight Capacity	Up to 350lbs
Seat Width	13" standard, 18" optional
Seat Height	20.5"
Battery	48V 12amp or optional 48V 6amp (travel)
Motor	48V -500W
Battery Charger	48V – 2AH