



QUALITY

SUPER SUPPORTIVE

Hand-blown, Overstuffed,
Conjugated Microfiber
Dacron Fill

EASY DISCONNECT

Quick-release Hand Control



STRONGEST IN THE INDUSTRY

- Northern Maple Hardwood Frame on all Stress Points
- Double Doweled at Stress Points
- Reinforced Mitre Blocks
- Stapled and Glued

LONG LASTING FOR INCREDIBLE COMFORT

5 Standard 8-gauge Springs

POCKETED COIL SEATING SYSTEM

Standard on most models*

(*See line sheet for applicable models)

FITS IN ANY HOME DÉCOR

Wide Fabric Selection

STABILITY & SUPPORT

Adjustable Anti-Skid Feet,
Rear Landing Legs**

(**Available on select models)

STRONG & LONG LASTING

9-gauge Steel
Scissor Mechanism



Complies with
UL 73 and
CSA C22.2 No. 68

SMALL TO TALL WE BUILD THEM ALL



Discover more at
www.goldentech.com