



# AL300RV

## RECREATIONAL VEHICLE LIFT



### REVOLUTIONARY BY DESIGN

Engineered to handle the stress of wind and heavy travel with robust hitch adapter to manage higher forces. Larger tubing in the platform avoids structural fatigue in off-road conditions.

### AL300RV PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb
ACCOMMODATES	<ul style="list-style-type: none"><li>- Power wheelchairs</li><li>- Scooters</li></ul>
HITCH CLASS	III or greater
POWERED	<ul style="list-style-type: none"><li>- Manual fold</li><li>- Power lifting</li></ul>
AVAILABLE OPTION(S)	<ul style="list-style-type: none"><li>- Battery pack</li><li>- Cover</li></ul>

INSTALLED WEIGHT	82 lb
PLATFORM SIZE	26" x 39" adjustable cradle extends up to 46"
WARRANTY	3-year
INCLUDED	<ul style="list-style-type: none"><li>- AL107 hitch post</li><li>- H300102 class III reinforced hitch adapter</li><li>- License plate mount</li></ul>
SECUREMENT	<ul style="list-style-type: none"><li>- 4 self-tensioning Q'Straint retractors</li><li>- Rear wheel cradle</li></ul>