

# mojo



## Mojo Auto-Fold

It's Freedom Baby, Yeah!

- 13 Mile Battery Range
- Airline Approved Battery
- Informative Digital Display
- Remote Operated Opening/Closing
- Lightweight & Durable Magnesium Alloy Construction
- Scooter Splits into Two Smaller Pieces for Even Lighter Lifting





The Mojo Auto-Fold conveniently opens and closes at the touch of a button.



Split the Mojo Auto-Fold in two for easier and lighter lifting.



The informative LCD Display shows the battery level, battery voltage, speed you are traveling at and distance.

**BUT WAIT!  
THERE'S MORE!**

The Mojo Auto-Fold features paddle acceleration, electromagnetic braking, LED headlight, basket, adjustable armrests, adjustable tiller, comfortable seating, and even a horn.



Enhance Your Life



Model	MOJO Auto-Fold (MJ101)
Overall Length	39 ¾"
Overall Width	18"
Tires	Honeycomb 6" Solid Front Tire
	Honeycomb 7" Solid Rear Tires
Maximum Speed	Up to 3.75 mph
Range	Up to 13 miles
Turning Radius	55"
Weight (without battery pack)	50 lbs.
Color	Red, Blue, Aqua, Lemony-Lime
Brake	Electromagnetic
Drive system	Rear Wheel Drive
Handlebar	PU Foam
Maximum weight capacity	Up to 300 lbs.
Seat Width	14"
Seat Height (from the ground)	23"
Battery	Dual (2) 24V 7.5AH
Motor	24V 120W
Battery Charger	29V 2AH
Folded Length	20 ½"
Folded Height	29"
Heaviest Piece	Front Frame 30.5 lbs.
Ground Clearance	1"

