

VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST COMFORTABLE RECLINER IN THE WORLD!

THE 4-MOTOR RECLINER



Shown in Badlands Mushroom Fabric



*Wireless Charging Station



COUNTLESS POSITION COMBINATIONS

ATLAS PLUS 2¹

When you need a moment for yourself, the **Atlas Plus 2** is there for you. Customize your comfort while you watch the big game or read your favorite book. Enjoy a refreshing beverage and keep it close with the cup holder, and conveniently charge your phone with the wireless charger.

The Atlas Plus 2 has been upgraded with enhanced true-infinite positioning and an ergonomically designed, memory toggle remote that features two programmable memory position buttons, as well as up/down/recline functions and customized comfort.



Scan here for more information on **VivaLift!®**









PRIDE®
MOBILITY
LIVE YOUR BEST®

ATLAS PLUS 2¹

VivaLIFT![®]
POWER RECLINER COLLECTION¹

COLLECTION	ATLAS PLUS 2 ¹
MODEL	PLR-2985M
POSITION	True Infinite
HEIGHT RANGE	5'4"-5'9"
SIZE	Medium
WEIGHT CAPACITY	400 lbs.
SEAT-TO-FLOOR	19.5"
SEAT DEPTH	21"
SEAT WIDTH	22"
TOP OF BACK TO SEAT	26"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	16"
OVERALL WIDTH	39"
CHAIR WEIGHT	164 lbs.
WARRANTY	7-Yr. Pro-rated

INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. **NEW!** Enhanced true-infinite position recline enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.

STANDARD FEATURES

- NEW!** Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions. Additionally, it has a convenient USB charger port that allows you to charge your smart device.
- Lithium battery backup ensures your recliner is operational when the power goes out.

1. VivaLift![®] Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

* The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.

STANDARD FABRIC: **BADLANDS**



MUSHROOM

WALNUT

STEEL

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]