



## Standard Colors:



Red

**ALL NEW!**  
**Great for Indoor &  
Outdoor Use!**

## Buzzaround 3 wheel GB107 4 wheel GB108



The all-new Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! Offered in both 3-wheel and 4-wheel models, this sporty design features a lightweight, stadium style seat, a weight capacity of 300 lbs., and increased foot room. The new Buzz LT also has a charger port in the tiller, a new hard plastic rugged basket, and color through shrouds. Enjoy its tight indoor maneuverability and smooth outdoor ride with 2.5" of ground clearance!

## Standard Features:

- Designed for indoor and outdoor use
- 300 lb. weight capacity
- Lightweight stadium style seat
- Increased foot room on a non-slip deck
- Durable color through shrouds
- Charging port on the tiller
- Disassembles quickly and easily for transport





## LIGHTWEIGHT PORTABLE SCOOTER

Model	GB107	GB108
Medicare Code (Pending)	K0800	K0800
Weight Capacity	300 lbs	300 lbs
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range <sup>1</sup>	7 miles	7 miles
Ground Clearance Under Deck	2.5"	2.5"
Turning Radius	39"	48"
Type Batteries	2-12AH	2-12AH
Colors	Red	Red
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	50 Amp	50 Amp
<b>Dimensions:</b>		
Length	41.75"	44"
Width	19-21"	19-21"
Height (ground to top of seat back)	32.5" - 34"	32.5" - 34"
Ground to Top of Seat	18.75" - 20.25"	18.75" - 20.25"
Ground to Top of Deck	4.25"	4.25"
Deck to Top of Seat	14.5" - 16"	14.5" - 16"
Number of Seat Height Positions	3	3
Size of Increments	.75"	.75"
Front Axle to Rear Axle	33.5"	33.5"
<b>Weight of Unit: (assembled)</b>		
Front Half	28 lbs	34 lbs
Rear Half	20.5 lbs	20.5 lbs
Seat w/Arms	13.5 lbs	13.5 lbs
Batteries (pack with 2 batteries)	17.5 lbs	17.5 lbs
<b>Tires:</b>		
	Solid, Low Profile	Solid, Low Profile
Front	8"	8"
Rear	8"	8"
Rear Anti-Tippers	2"	2"
<b>Standard Seat:</b>		
Back Height - no headrest	12.75"	12.75"
Width x Depth	16" x 14"	16" x 14"
Color	Black Vinyl	Black Vinyl
<b>Warranty:</b>		
Frame	Lifetime	Lifetime
Drive Train	2 Years	2 Years
Electronics <sup>2</sup>	2 Years	2 Years
Battery	1 Year	1 Year

<sup>1</sup> Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

<sup>2</sup> Electronics warranty excludes batteries.



Pivoting tiller folds easily.



Sleek dash with comfort grip handles and conveniently located charging port on the tiller.



Width-adjustable, flip up arm rests for customized comfort.



Quickly disassembles into five, lightweight pieces for easy transport!



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Find us on  
**Facebook**

[www.facebook.com/goldentechnologies](http://www.facebook.com/goldentechnologies)