

TRANSFORMER 2 FOLDING SCOOTER

Remote
Operated



- Contemporary New Look
- Folds Compact & Portable
- Up to 13.5 Miles Per Charge
- Lightweight
- Easy to Carry
- 300lb Weight Capacity
- Remote Operated
- Ideal for Travel
- Fits in Most Vehicles
- Airline Approved Battery
- Front & Rear Suspension
- LED Headlight



Shown in Blue
Available in Four Colors



Specifications	Transformer S3026
Overall Length	38"
Overall Width	18" w/ Armrests
Tires	Front – 6x1.5" Rear – 7x2.36"
Maximum Speed	3.75mph
Range	Up to 13.5 miles
Turning Circle	55"
Weight (w/ battery)	55lbs
Color	Red, Blue, Yellow, Black
Brake System	Intelligent Regenerative Electromagnetic Braking
Drive System	Rear Wheel Drive
Ground Clearance	3 1/4 "
Weight Capacity	Up to 300lbs
Seat Width	17"
Seat Height (from ground)	22"
Battery	24V10Ah Lithium Ion
Motor	120W
Battery Charger	DC24V/2A
Folded Dimensions	18" x 22" x 25"

TRANSFORMER 2