



## User & Vehicle Friendly

Ease of use and versatility makes the Bruno Curb-Sider® a time-tested favorite. Connect the docking device to a scooter or powerchair, push a button, and the Curb-Sider lifts, rotates and stores the mobility device. The Curb-Sider works with a wide variety of vehicles and mobility device weights.

**Automated.** Power lift/lower and rotate for effortless operation.

**Optimized.** Vehicle-specific bases provide precise installation. Telescoping head Curb-Sider model (VSL-6900) provides additional extension to clear bumpers and load wide mobility devices into smaller openings.

**Slim.** Preserve maximum seating and storage. Partial third-row seating or cargo space is available when a mobility device is not present. Optional quick-release pins enable the lift head to fold down for even greater space.

**Adaptable.** Reinstall in another applicable vehicle or adjust to accommodate a new mobility device.

### Lifting Capacity

250 lb (114 kg)

400 lb (181 kg)

450 lb (203 kg)

### Vehicle Fit

SUVs, Vans, Pickups

### Warranty

3-Year Limited

### Learn More

[Bruno.com](http://Bruno.com)



## OPTIONS

- Exterior control VSL-6000 only. (A)  
(Included with VSL-6000 450lb)
- Quick release pins to fold down lift head. (B)
- Driver-side placement (C)
- Back-Off Device easily removes back of mobility device. (D)

