

Independence Freedom Mobility

zeen[®]
by exokinetics

The Zeen is an innovative mobility device filling the gap between a walker and a wheelchair. Providing a power-free lift between sitting and standing, you can enjoy effortless walking with the safety and reassurance of support through your hips and buttocks at all times.

The Zeen is for the person who values independent mobility in a standing position and is experiencing trouble with strength, balance, endurance or fear of falling.

**Let us safely get you
on your feet, stay there
with confidence and
live your best life!**



Is the Zeen right for me?

- Do you have a strong desire to be independent with mobility?
- Are you able to walk multiple steps with or without assistance?
- Do you have moderate grip strength in at least one hand?
- Is rising from a seated position challenging?
- Have you limited walking/standing due to fear of falling, strength, balance or endurance challenges?
- Are you able to learn to use a new device on your own, or do you have someone to assist you?

If you answered YES to any of these questions, the Zeen could be the mobility device for you!

HEALTH BENEFITS ASSOCIATED WITH STANDING AND EYE-LEVEL INTERACTION

- Improves bone & joint integrity, circulation, leg strength
- Strengthens cardiovascular and digestive systems
- Decreases risk of heart disease, diabetes, cancer, low back pain, weight gain and obesity
- Improves well-being, alertness, mood, thinking, cognition, sleep quality
- Increases energy level, balance and coordination
- Decreases stress, fatigue, depression and anxiety

 **sit. stand. go**

Zeen by Exokinetics
19 Hagerty Blvd. Unit 15, West Chester, PA 19382



gozeen.com




“It has been amazing — I am able to walk and stand up straight without the fear of falling. I know that whenever I needed a few minutes to rest, I have a seat right behind me.” — Amy G.

Co-inventor Garrett Brown, renowned in the filmmaking industry for his creation of the Steadicam, watched his parents and their friends grapple with declining mobility. He envisioned a mobility device that would allow people to “float” to standing and move easily through any room, hallway, or outdoor setting under their own power.

ZEEN MINIMUM USER REQUIREMENTS

- Powerful desire to be independent and mobile
- Able to stand and walk short distances with or without assist
- Ability to learn and retain new information
- Weight between 50-250 lbs. (23-113 kg)
- Distance from ground to groin between 21-35" (53-89cm)
- Moderate grip strength in at least 1 hand

PRODUCT NAME	SKU
Zeen Standard	
- Foam Wheels	ZEEN-000
- Pneumatic Wheels	ZEEN-002
Zeen One-Hand Control	
- Foam Wheels	ZEEN-004
- Pneumatic Wheels	ZEEN-002
Spare Wheel Kit	
- 2 Never Flat Wheels	100002
- 2 Pneumatic Wheels	100001
Cup Holder	002-A9004-50
Travel Cover	002-A9005-60
Trunk Lift Straps	002-A9002-60
Travel / Exercise Pins	1000003



Zeem folds for easy transport



SPECIFICATIONS

Dimensions (LxWxH)	FOLDED 33.5" x 11" x 38"
SEATED MODE 33.5" x 25.2" x 36.8"	LOWEST SEAT HEIGHT 20.5"
MAX STANDING HEIGHT 33.5" x 25.2" x 51.3"	WEIGHT ~42 pounds

CLEARANCE

Doorways	27"
Turning diameter	45"
Width of seat between armrests	21"

DEALER INFO



Covering Florida

P: (813) 642-0322

E: info@mymobilitystore.com

**729 Cypress Village Blvd,
Ruskin, FL 33573**

www.mymobilitystore.com