




Standard Fabrics:




Calypso



Copper



Hazelnut



Shiraz



Sterling

This fabric stands out for its comfort, smooth texture, and rich color options. It features our exclusive SofTouch™ with Stain Defense coating.

* An upcharge will apply for all crossover fabrics.

Cambridge PR401



Comfort Zone[™]
by GOLDEN

The popular Cambridge power lift recliner features our hybrid sizing for a better fit. The Small/Medium size fits people 5'1" to 5'6" and the Medium/Large size fits people 5'7" to 6'2". The Cambridge features a supportive biscuit backrest design and a soft, comfortable seat for hours of relaxation.

Standard Features:

- Classic Design
- Blissful Comfort with Ultra Soft Back
- Simple Two-Button Hand Control
- USB Charging Port



SIGNATURE Series PR401

Model Number	PR401-SME	PR401-MLA
Size	Small/Med	Med/Large
Weight Capacity	375 Lb.	375 Lb.
Positioning	3	3
Comfort Zone(s)	1	1
Twilight Technology	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	NO	NO
Back Type	Biscuit	Biscuit
Removable Back Standard	YES	YES
Head Rest Cover	YES	YES
Arm Rest Covers	YES	YES
Overall Width	33"	35"
Overall Height	41.5"	44.5"
Width Between Arms	20"	22"
Floor to Top of Seat	19"	22"
Seat Depth	19"	22"
Seat to Top of Back	27"	31"
Seat to Top of Arm	6.5"	7"
Distance Required from Wall Reclined	18"	19"
Extension reclined from Back of Seat to End of Foot Rest	41.5"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Pocket	Inside Pocket
Extra Pocket Option	YES	YES
Heat & Massage Option	YES	YES
Articulating Headrest Option	NO	NO
Adjustable Lumbar Option	NO	NO
AutoFlex Option	NO	NO
Foot Rest Extension Option	YES	YES
Straight Lift Option	YES	YES
USB Charging Port	YES	YES
Battery Backup	YES	YES
Place of Origin	USA	USA



Upgrade with
Deluxe Therapeutic
Heat & Massage!



Shown in Calypso


SmartTek Diagnostic System STANDARD on all Golden models.

Comfort Zone™ by GOLDEN

Comfort Zone 1 provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.



WE PATENTED COMFORT™



Scan with your
smartphone
camera to learn
more about
Comfort Zone
by Golden.



Phone: 800-624-6374
Fax: 800-628-5165
www.goldentech.com



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



Find us on
Facebook
www.facebook.com/goldentechnologies

SHT-PR401-111021