




Shown in Maple

Standard Fabrics:



Acorn



Graphite



Maple

This faux leather is a rich and modern alternative to natural leather that is extremely comfortable while offering the rugged durability of leather. Made from 100% polyester, the fabric is easy to care for ensuring functionality as well as elegance. This long-lasting fabric will complement a variety of decors while standing the test of time under heavy use.

Crossover fabrics are unavailable.

Dione PR446



Comfort Zone[™]
by GOLDEN

The Dione is offered in both medium and large sizes and features a stylish, overstuffed seamed back and plush chaise pad. It comes standard with adjustable lumbar and headrest, programmable AutoDrive[™] hand control, and USB charging port for your phone or tablet. Also standard is our patented four Comfort Zone positioning system, which provides a stress-free and supportive reclining experience.

Standard Features:

- Patented MaxiComfort[®] Positioning
- Adjustable Headrest
- Adjustable Lumbar
- Programmable AutoDrive[™] Hand Control
- USB Charging Port
- Extra Pocket



DeLUNA Series PR446

| | | | | | | | | | | | | | | | | | | | | |
|--|-------------|--|-----------------|----|----|----------------------------|----|----|----------------------|-----|-----|-------------------|-------------|-------------|----------------|----|----|-----------------|-------|-------|
| Model Number | PR446-MED | PR446-LAR | | | | | | | | | | | | | | | | | | |
| Size | Medium | Large | | | | | | | | | | | | | | | | | | |
| Weight Capacity | 375 LBS | 375 LBS | | | | | | | | | | | | | | | | | | |
| Positioning | MaxiComfort | MaxiComfort | | | | | | | | | | | | | | | | | | |
| Comfort Zone(s) | 4 | 4 | | | | | | | | | | | | | | | | | | |
| Twilight Technology | NO | NO | | | | | | | | | | | | | | | | | | |
| Chaise Seat | YES | YES | | | | | | | | | | | | | | | | | | |
| Pocketed Coil Seat | YES | YES | | | | | | | | | | | | | | | | | | |
| Back Type | Seam | Seam | | | | | | | | | | | | | | | | | | |
| Removable Back Type | Bracket | Bracket | | | | | | | | | | | | | | | | | | |
| Head Rest Cover | NO | NO | | | | | | | | | | | | | | | | | | |
| Arm Rest Covers | NO | NO | | | | | | | | | | | | | | | | | | |
| Overall Width | 34" | 35" | | | | | | | | | | | | | | | | | | |
| Overall Height | 43.5" | 46" | | | | | | | | | | | | | | | | | | |
| Width Between Arms | 22" | 22.5" | | | | | | | | | | | | | | | | | | |
| Floor to Top of Seat | 20.5" | 22" | | | | | | | | | | | | | | | | | | |
| Seat Depth | 21" | 21.5" | | | | | | | | | | | | | | | | | | |
| Seat to Top of Back | 28" | 30" | | | | | | | | | | | | | | | | | | |
| Seat to Top of Arm | 7" | 7.5" | | | | | | | | | | | | | | | | | | |
| Distance Required from Wall Reclined | 24" | 26" | | | | | | | | | | | | | | | | | | |
| Extension Reclined from Back of Seat to End of Foot Rest | 42.5" | 44" | | | | | | | | | | | | | | | | | | |
| Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism | YES | YES | | | | | | | | | | | | | | | | | | |
| Warranty on Electrical Parts | 3 YR | 3 YR | | | | | | | | | | | | | | | | | | |
| Pro Rated Years 4 Thru 7 | YES | YES | | | | | | | | | | | | | | | | | | |
| Warranty on Mechanical Labor | 3 YR | 3 YR | | | | | | | | | | | | | | | | | | |
| Warranty on Electrical Labor | 1 YR | 1 YR | | | | | | | | | | | | | | | | | | |
| Hand Control Routing | Inside Arm | Inside Arm | | | | | | | | | | | | | | | | | | |
| Extra Pocket Option | Standard | Standard | | | | | | | | | | | | | | | | | | |
| Heat & Massage Option | NO | NO | | | | | | | | | | | | | | | | | | |
| Articulating Headrest Option | Standard | Standard | | | | | | | | | | | | | | | | | | |
| Adjustable Lumbar Option | Standard | Standard </tr <tr> <td>AutoFlex Option</td> <td>NO</td> <td>NO</td> </tr> <tr> <td>Foot Rest Extension Option</td> <td>NO</td> <td>NO</td> </tr> <tr> <td>Straight Lift Option</td> <td>YES</td> <td>YES</td> </tr> <tr> <td>USB Charging Port</td> <td>Outside Arm</td> <td>Outside Arm</td> </tr> <tr> <td>Battery Backup</td> <td>NO</td> <td>NO</td> </tr> <tr> <td>Place of Origin</td> <td>China</td> <td>China</td> </tr> | AutoFlex Option | NO | NO | Foot Rest Extension Option | NO | NO | Straight Lift Option | YES | YES | USB Charging Port | Outside Arm | Outside Arm | Battery Backup | NO | NO | Place of Origin | China | China |
| AutoFlex Option | NO | NO | | | | | | | | | | | | | | | | | | |
| Foot Rest Extension Option | NO | NO | | | | | | | | | | | | | | | | | | |
| Straight Lift Option | YES | YES | | | | | | | | | | | | | | | | | | |
| USB Charging Port | Outside Arm | Outside Arm | | | | | | | | | | | | | | | | | | |
| Battery Backup | NO | NO | | | | | | | | | | | | | | | | | | |
| Place of Origin | China | China | | | | | | | | | | | | | | | | | | |



Shown in Graphite



Convenient built-in USB charging port for your phone or tablet located on the outside arm.

SmartTek Diagnostic System STANDARD on all Golden models.
Designed and engineered in Old Forge, PA. Assembled in China using Golden components.




Comfort Zone™ by GOLDEN

Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374
Fax: 800-628-5165
www.goldentech.com

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.

