

# Lightweight mobility solutions

Rediscover your freedom





# Recognised leaders in the design and manufacture of ultra-portable scooters and power chairs

MK1.5

Lite

Power Chair

**My name is Sumi, and welcome to eFOLDi - Home to the world's lightest and most compact folding mobility scooter and power chair.**

We strive to help people rediscover their freedom and independence.

When my father lost his mobility in a work-place accident some years ago, he simply couldn't find a mobility scooter anywhere, that was easy to use, safe and comfortable - but also lightweight and portable. Undeterred, my father set about designing his own scooter – and that is how the first eFOLDi mobility scooter came to be born!

Since then, eFOLDi has developed into a thriving family business, designing and manufacturing innovative mobility products that help make travel simple, achievable and most importantly, fun! The eFOLDi scooter was also a proud winner of the prestigious 'British Invention Of The Year' award in 2016. Legendary innovator and entrepreneur, Sir Richard Branson, even named the eFOLDi scooter "folding magic".

Join us in making travel easy and enjoyable again, and let an eFOLDi be your perfect travel companion.

As seen on TV

*Sumi Wang*  
**Founder**



**FREEDOM • FRIENDSHIP • FUN**

## What makes eFOLDi unique



### Lightweight

Carry and lift with ease



### Portable

Easy to take with you, anywhere



### Compact

Folds down to cabin luggage size.



### Robust

Peace of mind with strong, durable materials.



*"eFOLDi is folding magic."*

*Sir Richard Branson*



MK1.5



## Ultra-lightweight folding scooter

The multi award winning MK1.5 from eFOLDi is the world's lightest, most compact folding UK road-legal scooter. This ultra-portable scooter folds and unfolds simply, in a matter of seconds and requires no complex assembly - just unfold and go.

Lightweight • Portable • Compact • Robust



**Lightweight**  
at only 17kg/37.5lbs\*



**Easy to use**  
controls



**Compact size**  
folds with ease



**Airline friendly battery**  
14 miles/22km  
per charge



**Travel in comfort**  
suspension/  
pneumatic tyres



**Fully UK road legal**  
8mph/12kmh  
top speed

“

*My eFOLDi has impacted hugely on my life. It has enabled me to go out so much more than before. Now I have got my freedom and independence back, I can go out whenever and wherever I want.*

**Ray, London**

”

\*excl. battery

Lite



## Ultra-lightweight folding scooter

The eFOLDi Lite is the world's lightest, most compact folding travel-scooter. The Lite is the latest, and most advanced scooter from eFOLDi - weighing just 15kg\*. It folds and unfolds in just a few seconds and is easy to store and transport. No complex assembly is required - just unfold and go.

**Lightweight • Portable • Compact • Robust**



**Lightweight**  
at only 15kg/33lbs\*



**Easy to use**  
thumb controls



**Compact size**  
folds with ease



**Travel in comfort**  
suspension/  
pneumatic tyres



**Airline friendly battery**  
10 miles/16km  
per charge



**Robust**  
and carries up to  
19 Stone/120kg/265lbs

\*excl. battery



*eFOLDi has changed my life so much for the better. It allows me to participate so much more in the world. It folds up extremely neatly like a small suitcase, and is light enough for someone to easily lift into a car boot. I have gone from being largely isolated to being much more able to participate in outings.*

**Jackie, Esher**



## Power Chair



## Ultra-lightweight folding power chair

The eFOLDi Power Chair is the world's lightest, most compact folding powered wheelchair. This ultra-portable, lightweight power chair is the perfect travel companion, and folds and unfolds in a matter of seconds. No complex assembly is required - just unfold and go.

**Lightweight • Portable • Compact • Robust**





**Handy controls**  
for attendant to drive



**Easy to use**  
joystick and control pad



**Compact size**  
folds with ease



**Travel in comfort**  
durable, breathable fabric



**Airline friendly batteries**  
8 miles/13km per charge



**Lightweight**  
at only 15kg/33lbs\*

\*excl. battery



*With an eFOLDi you can just fold it up, and take it with you on buses, trains, or planes. Perfect for me to get out and about with my family and even go on holiday abroad!*

**Anna, Reigate**



# Specifications

MK1.5



Lite



Power Chair



Maximum Speed	8mph / 12kph	4mph / 6kph	4mph / 6kph
Maximum Range	14 miles / 22km	10 miles / 16km	8 miles / 13km
Maximum User Weight	19 Stone / 120kg / 265lbs	19 Stone / 120kg / 265lbs	21 Stone / 135kg / 297lbs
Climbing Ability	1:6 Gradient	1:10 Gradient	1:10 Gradient
Folded Dimensions	26 x 22 x 16 inches 66 x 56 x 40cm	26 x 22 x 15 inches 65 x 55 x 37cm	22 x 13 x 36 inches 58 x 34 x 92cm
Weight (excl. battery)	37.5 lbs /17kg	33lbs /15kg	33lbs /15kg
Weight (incl. battery)	42lbs/19kg	37.5lbs/17kg	37.5lbs/17kg
Type of Use	Road (UK legal)/Pavement	Pavement Only	Pavement Only
Braking System	Manual	Electro-Magnetic	Electro-Magnetic
Airline friendly	Yes	Yes	Yes

*Specifications may vary due to manufacturing tolerances.*

## Accessories

To complement your eFOLDi, we offer a wide range of accessories tailored to suit a variety of needs. This includes a host of storage solutions, protective covers, cushions and battery charger options.



## Innovation

**Recognised as leaders in the design and manufacture of ultra-portable scooters and power chairs, we strive to help people rediscover their freedom and independence.**

Our products are manufactured using the highest quality materials to ensure they're robust enough to withstand the toughest of conditions and terrain. Our scooter frames are made from high-grade aluminium or magnesium alloys, our high-specification brushless motors deliver exceptional

performance, and our batteries use advanced lithium-ion technology, so they're not only durable, they offer outstanding levels of reliability.



## Support

All eFOLDi products come with a comprehensive 12-month UK manufacturer's warranty. In the unlikely event of a problem, our UK support centre and servicing team are always on hand to help. They are trained to diagnose and rectify any problems, including warranty repairs and any ongoing service and maintenance requirements.



**For more information...**

**Call us today on:  
0800 999 8123**

**Or email us at:  
sales@efoldi.com**



Authorised eFOLDi Dealer:

eFOLDi 

0800 999 8123 • [sales@efoldi.com](mailto:sales@efoldi.com) • [www.efoldi.com](http://www.efoldi.com)

SunTech UK Limited t/a eFOLDi, 25 Ormside Way, Holmethorpe Industrial Estate, Redhill, Surrey RH1 2LW.



**FREEDOM • FRIENDSHIP • FUN**